

Nitinol Tipless Baskets

HALO™ 1.5 Fr Sheath Diameter

- Exceptional irrigation enhances visibility.
- Enhanced kink resistance increases longevity.
- Can share 3.6 Fr channel with a laser fiber.

Vantage™ 2.4 Fr Sheath Diameter





- Maintains basket symmetry.
- Enhances capture of small fragments.
- Minimizes tissue trauma.
- Increased rigidity for larger stone fragments.

Halo Trinity™ 1.5 Fr Sheath Diameter

- Wide proximal gaps for easy stone engagement.
- Able to extract multiple stone fragments simultaneously.

6-Wire HALO™ 1.5 Fr Sheath Diameter


- Effectively captures both small and large stones.
- Able to extract multiple stone fragments simultaneously.



Helical Baskets

HERCULES™ Stainless Steel Paired Helical Basket


- Maximum dilation force.
- Available with or without filiform tip.



Flat Wire Baskets

APEX™ Nitinol Flat Wire Basket

- Maximum basket durability.



Handles

- Allows independent basket rotation.
- Ultra lightweight and highly tactile.



Nitinol Tipless Baskets

ORDER NO.	SHEATH DIA. (Fr)	SHEATH LENGTH (cm)	SHEATH MATERIAL	BASKET DIA. (mm)	NO. OF WIRES	TIP	COMMENTS
HALO™ Nitinol Tipless Basket							
554120	1.5	120	Polyimide	12	4	Tipless	Optimal Irrigation
6-Wire HALO™ Nitinol Tipless Basket							
556160	1.5	120	Polyimide	16	6	Tipless	
Halo Trinity™ Nitinol Tipless Basket							
556130	1.5	120	Polyimide	11	6	Tipless	
Vantage™ Nitinol Tipless Basket							
550320	2.4	120	Polyimide	16	4	Tipless	

Helical Baskets

ORDER NO.	SHEATH DIA. (Fr)	SHEATH LENGTH (cm)	SHEATH MATERIAL	BASKET DIA. (mm)	NO. OF WIRES	TIP	COMMENTS
HERCULES™ Stainless Steel Paired Helical Basket							
530520	3.0	120	Polyimide	11	2 Singles, 2 Pair	"O" Tip	Maximum Dilation Force
530525	3.0	120	Polyimide	11	2 Singles, 2 Pair	5cm filiform	

Flat Wire Baskets

ORDER NO.	SHEATH DIA. (Fr)	SHEATH LENGTH (cm)	SHEATH MATERIAL	BASKET DIA. (mm)	NO. OF WIRES	TIP	COMMENTS
APEX™ Nitinol Flat Wire Basket							
544310	2.4	90	Polyimide	15	4	"O" Tip	
544320	2.4	120	Polyimide	15	4	"O" Tip	

Caution: Federal (USA) law restricts these devices to sale by or on the order of a physician.

Sacred Heart Medical, Inc.

Distributed By:

6018 Blue Circle Drive
Minnetonka, MN 55343 USA

Tel: (952) 345-0280
Fax: (952) 942-7590
contact@sacredheartmedical.com
www.sacredheartmedical.com

©2025 Sacred Heart Medical, Inc.
All rights reserved. 93-7001-01 Rev. A (2025-02)