



Stone Baskets

Nitinol Tipless Baskets

- Baskets designed to maintain symmetry
- Kink resistant nitinol wire


HALO™ 1.3 or 1.5 Fr Sheath Diameter

- Small sheath diameter allows exceptional irrigation flow to enhance visibility
- Can share 3.6 Fr working channel with a laser fiber




Halo Trinity™ 1.5 Fr Sheath Diameter

- Wide proximal gaps for stone engagement
- 6 wires on distal end to retain small fragments



6-Wire HALO™ 1.5 Fr Sheath Diameter

- 6 wire design for secure stone capture
- 16 mm basket diameter





Vantage™ 2.4 Fr Sheath Diameter

- Increased rigidity and basket dilation force
- 16 mm basket diameter

Helical Baskets

HERCULES™ 3.0 Fr Sheath Diameter


- Maximum basket dilation force
- Available with or without filiform tip



Flat Wire Baskets

APEX™ 2.4 Fr Sheath Diameter

- Kink resistant nitinol wire



Ergonomic Handle

- Independent basket rotation feature
- Ultra lightweight and highly tactile



Nitinol Tipless Baskets

ORDER NO.	SHEATH DIA. (Fr)	SHEATH LENGTH (cm)	BASKET WIRE MATERIAL	BASKET DIA. (mm)	NO. OF WIRES	TIP
HALO™ Tipless Stone Basket						
554120	1.5	120	Nitinol	12	4	Tipless
555020	1.3	120	Nitinol	12	4	Tipless
Halo Trinity™ Tipless Stone Basket						
556130	1.5	120	Nitinol	11	6	Tipless
6-Wire HALO™ Tipless Stone Basket						
556160	1.5	120	Nitinol	16	6	Tipless
Vantage™ Tipless Stone Basket						
550320	2.4	120	Nitinol	16	4	Tipless

Helical Baskets

ORDER NO.	SHEATH DIA. (Fr)	SHEATH LENGTH (cm)	BASKET WIRE MATERIAL	BASKET DIA. (mm)	NO. OF WIRES	TIP
HERCULES™ Helical Stone Basket						
530520	3.0	120	Stainless Steel	11	2 Singles, 2 Pair	"O" Tip
530525	3.0	120	Stainless Steel	11	2 Singles, 2 Pair	5 cm filiform

Flat Wire Baskets

ORDER NO.	SHEATH DIA. (Fr)	SHEATH LENGTH (cm)	BASKET WIRE MATERIAL	BASKET DIA. (mm)	NO. OF WIRES	TIP
APEX™ Flat Wire Stone Basket						
544310	2.4	90	Nitinol	15	4	"O" Tip
544320	2.4	120	Nitinol	15	4	"O" Tip

Intended Purpose: The Stone Basket is intended for endoscopic manipulation and removal of stones in the urinary tract.

Some products may not be available in all markets

Caution: Federal (USA) law restricts these devices to sale by or on the order of a physician.

Sacred Heart Medical, Inc.

Distributed By:

6018 Blue Circle Drive
Minnetonka, MN 55343 USA

Tel: (952) 345-0280
Fax: (952) 942-7590
contact@sacredheartmedical.com
www.sacredheartmedical.com

©2026 Sacred Heart Medical, Inc.
All rights reserved. 93-7001-01 Rev. B (2026-06)